

Simultaneous HIS and
RVA Recording with a
Single Catheter

APT-EP Map-iT[®] HIS-RV[™] Diagnostic Catheters



KEY FEATURES



- > Decapolar catheters with four distal electrodes for the RVA and six proximal for the HIS
- > Simultaneous HIS and RVA recording with a single catheter
- > Pace the RVA and record the retrograde AV node activation with a single catheter
- > Achieve simplified, more cost-effective procedures
- > With the superior shaft construction, maintains stable curve and shaft performance

CATHETER SPECIFICATIONS

- > **French Size:** 5F or 6F
- > **Length:** 120cm
- > **Electrodes:** 10: four distal (RVA) electrodes with six proximal (HIS area)
- > **Electrode Spacing:**
Variable gap between HIS and RV:
20, 23, 24, or 29mm
- > **Curve Type:** MPA



APT EP MAP-IT® HIS-RV™ DIAGNOSTIC CATHETERS

APT EP Map-iT® 5F Decapolar Fixed Curve Electrophysiology Catheter

Reference Number	French	Electrodes	Electrode Spacing (mm)	Ring Electrode (mm)	Tip	Curve	Useable Length (cm)
901578	5Fr	10	4-4-4-29-3-3-3-3-3	1	2	MPA	120
901579	5Fr	10	4-4-4-23-4-4-4-4-4	1	2	MPA	120
901580	5Fr	10	4-4-4-29-4-4-4-4-4	1	2	MPA	120
901581	5Fr	10	4-4-4-23-5-5-5-5-5	1	2	MPA	120
901582	5Fr	10	4-4-4-29-5-5-5-5-5	1	2	MPA	120
901583	5Fr	10	5-5-5-20-3-3-3-3-3	1	2	MPA	120
901584	5Fr	10	5-5-5-24-3-3-3-3-3	1	2	MPA	120
901585	5Fr	10	5-5-5-20-4-4-4-4-4	1	2	MPA	120
901586	5Fr	10	5-5-5-24-4-4-4-4-4	1	2	MPA	120

Connection Cable: Reference number 901957



APT EP Map-iT® 6F Decapolar Fixed Curve Electrophysiology Catheter

Reference Number	French	Electrodes	Electrode Spacing (mm)	Ring Electrode (mm)	Tip	Curve	Useable Length (cm)
901822	6Fr	10	4-4-4-29-3-3-3-3-3	1	2	MPA	120
901823	6Fr	10	4-4-4-23-4-4-4-4-4	1	2	MPA	120
901824	6Fr	10	4-4-4-29-4-4-4-4-4	1	2	MPA	120
901825	6Fr	10	4-4-4-23-5-5-5-5-5	1	2	MPA	120
901826	6Fr	10	4-4-4-29-5-5-5-5-5	1	2	MPA	120
901827	6Fr	10	5-5-5-20-3-3-3-3-3	1	2	MPA	120
901828	6Fr	10	5-5-5-24-3-3-3-3-3	1	2	MPA	120
901829	6Fr	10	5-5-5-20-4-4-4-4-4	1	2	MPA	120
901830	6Fr	10	5-5-5-24-4-4-4-4-4	1	2	MPA	120

Connection Cable: Reference number 901957