

HIS RV™ DIAGNOSTIC CATHETER

Multi-Electrode Catheter for Simultaneous HIS & RV Recording



Intentional Electrode Placement

Receive high quality electrogram data from four distal RV and six proximal HIS electrodes.



Streamlined Procedure

With a single catheter, pace the RV and record the AV node activation for a more efficient procedure.



Enhanced Access & Patient Safety

Reduce load on vasculature and preserve the opportunity for additional catheter access.


The Map-iT HIS RV, also known as the HIS-RV Diagnostic Catheter, allows physicians to use one catheter and a single access point for both HIS and RV recordings, simplifying procedures and potentially reducing the risk of HV interval prolongation. Map-iT catheters are designed with a braid-reinforced shaft for high pushability, easy placement, and consistent handling, even during long procedures.

Specifications

French Size	5F or 6F
Length	120cm
Electrodes	4 Distal, 6 Proximal
HIS to RV Electrode Spacing	20, 23, 24, 29mm
Curve Type	MPA



Superior Shaft Construction
High-performance torque control, pushability, and durability



5F Decapolar Fixed Curve Catheter

French	5F
Electrodes	10
Ring Electrode	1mm
Tip Electrode	2mm
Curve	MPA
Usable Length	120cm

Reference No.	Electrode Spacing (mm)
901579	4-4-4-23-4-4-4-4-4
901580	4-4-4-29-4-4-4-4-4
901583	5-5-5-20-3-3-3-3-3
901584	5-5-5-24-3-3-3-3-3

Connection Cable Reference No. 901957

Made in the USA

6F Decapolar Fixed Curve Catheter

French	6F
Electrodes	10
Ring Electrode	1mm
Tip Electrode	2mm
Curve	MPA
Usable Length	120cm

Reference No.	Electrode Spacing (mm)
901823	4-4-4-23-4-4-4-4-4
901824	4-4-4-29-4-4-4-4-4
901827	5-5-5-20-3-3-3-3-3
901828	5-5-5-24-3-3-3-3-3

Connection Cable Reference No. 901957

Made in the USA



To Order

✉ Map-iT@stereotaxis.com ☎ +1 (763) 432-9004

Follow Us

✉ in EP Map-iT